

ABRASION & PILLING TESTER

MARTINDALE



For the control of abrasion and pilling on almost any kind of woven and knitted fabrics, non-wovens, socks, gloves, natural and artificial leather, both dry and wet samples.

Extremely versatile, robust, accurate, silent, vibration free, 100% made in Italy.

Supplied with a basic starter kit of top quality consumables complete with traceable certificate.

ISO 17025 Calibration Certificate (Accredia - ILAC) available on demand.



COMPANY WITH
MANAGEMENT SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =
= ISO 14001 =

MARTINDALE

Code 2568 (9 positions) - Code 2568A (6 positions)



Compulsory Accessories for 6- & 9-position setting (at least 1 is needed)

Set of 6 or 9 sample holders Ø 38 mm to perform ABRASION test according to UNI EN ISO 12947-1, and to perform PILLING test according to ASTM D4970, Code 2568A.900, or Code 2568.900.

Set of 6 or 9 sample holders Ø 90 mm to perform PILLING test according to UNI EN ISO 12945-2, Code 2568A.300, or Code 2568.300.

Optional Accessories (available on demand)

Sample holder for socks & hosiery ABRASION test according to UNI EN 13770, Code 2568.804.

Sample holder for ABRASION test on leather and on coated fabrics, according to UNI EN ISO 17076-2, Code 2568.260.



Example of 1-position sample holder with Ø 38 mm



Example of 1-position sample holder with Ø 90 mm



Example of 1-position sample holder for socks & hosiery



Example of 1-position sample holder for leather and coated fabrics



Testing procedure

In conformity with the International Standards, the fabric specimen to be tested is brought into contact and rubbed against a Standard reference fabric (as regards abrasion test), or against the same type of specimen, or against a standard felt (as regards Pilling test), following a predetermined number of specific rotary motions known as “Lissajous figure” paths.

Basically, the abrasion test end point is reached once the operator visually detects the first two broken threads on the tested fabric specimen, whereas the Pilling test is accomplished by observing (after a certain number of laps), the created pill effect (fuzzing surface / surface disturbance) on the fabric specimen and comparing it with the special photographic standards.

For Pilling tests according to EN, ISO, and EMPA (also called “Swiss Pilling”) Standards, a particular sample holder is required (with special dimensions and weight) as well as special photographic standards: EMPA 991 for woven fabrics, and EMPA 992 for knitted fabrics.

Technical Features

Two models available: Martindale 6 positions, Code 2568A; Martindale 9 positions, Code 2568.

Modular & Flexible: the 6-position model can be further upgraded up to 9 positions by means of specific upgrading kits.


User-friendly control panel characterised by:

- wide LCD touch screen display for easy setting of testing parameters;
- graphic layout of the testing heads;
- graphic scheme to select the type of test and the movement path;
- start/stop and emergency keys.




Three types of test movement:

abrasion (big Lissajous), pilling (small Lissajous) and rubbing (linear). The desired movement is performed according to the set position of the 3 driving pivots, as shown by the graphic scheme on the control panel.



Intuitive software for the setting of testing parameters, such as number of cycles per testing step (with consequent stop), speed (23,75 - 47,5 - 71,25 rpm), increase/decrease counting modality, display language.



Easy access to all the testing positions.

Sample and abrasion surface can be easily inspected with no need to remove the movable plate.

ISO 9001 Calibration Report available on request.

ISO 17025 Calibration Certificate (Accredia - ILAC) available on demand.



MARTINDALE

AVAILABLE MODELS

Martindale 6 positions	Code	2568A
Martindale 9 positions	Code	2568

COMPULSORY ACCESSORIES FOR 6-POSITION SETTING (at least one is needed)

Set of 6 sample holders Ø 38 mm for ABRASION test according to UNI EN ISO 12947-1 and PILLING test according to ASTM D4970	Code	2568A.900
Set of 6 sample holders Ø 90 mm for PILLING test according to UNI EN ISO 12945-2	Code	2568A.300

COMPULSORY ACCESSORIES FOR 9-POSITION SETTING (at least one is needed)

Set of 9 sample holders Ø 38 mm for ABRASION test according to UNI EN ISO 12947-1 and PILLING test according to ASTM D4970	Code	2568.900
Set of 9 sample holders Ø 90 mm for PILLING test according to UNI EN ISO 12945-2	Code	2568.300

COMMON OPTIONAL CONSUMABLES FOR MARTINDALE 6 & 9 POSITIONS

1 roll of standard abrasive cloth, 160x100 cm for ABRASION test according to UNI EN ISO 12947-1	Code	314.12
24 discs Ø 140 mm of standard woven pre-cut felt for ABRASION test according to UNI EN ISO 12947-1	Code	314.8
25 pcs of Standard polyurethane foam, 25x20 cm for ABRASION test according to UNI EN ISO 12947-1	Code	314.32
24 discs Ø 90 mm of standard woven felt for PILLING test according to UNI EN ISO 12945-2	Code	314.20
24 discs Ø 140 mm of standard backing non-woven pre-cut felt for ABRASION test according to UNI EN ISO 12947-1-2-3-4 and M&S Standards	Code	314.26
Set of abrasive paper (25 pcs.) for ABRASION test according to UNI EN 388 (protective clothing)	Code	2560.326
Bi-adhesive tape for ABRASION test according to UNI EN 388 (protective clothing)	Code	2560.328
Set of abrasive paper (25 pcs.) for ABRASION test according to UNI EN ISO 5470-2 (coated fabric)	Code	2560.330
Cork base for sample cutters	Code	175B.46

COMMON STANDARD PHOTOGRAPHS FOR 6 & 9 POSITIONS - OPTIONAL

Set of comparative photographs for PILLING test of woven fabrics according to UNI EN ISO 12945 (EMPA 991)	Code	314.18
Set of comparative photographs for PILLING test of knitted fabrics according to UNI EN ISO 12945 (EMPA 992)	Code	314.24
Set of comparative photographs SM50 for PILLING test of woven fabrics according to ASTM D4970 and IWS TM 196	Code	314.14
Set of comparative photographs SM54 for PILLING test of knitted fabrics according to ASTM D4970 and IWS TM 196	Code	314.16

COMMON SAMPLE CUTTERS FOR 6 & 9 POSITIONS - OPTIONAL

Sample cutter to obtain Ø 140 mm discs for PILLING test	Code	2560.320
Sample cutter to obtain Ø 38 mm discs for ABRASION & PILLING test	Code	2560.322
Sample cutter to obtain Ø 90 mm discs for PILLING test, if pre-cut discs are not available	Code	2560.324

COMMON OPTIONAL FOR 6 & 9 POSITIONS - OTHER ACCESSORIES

Writing pivot to check the Lissajous figure	Code	2568.222
Sample holder for ABRASION test on socks & hosiery according to UNI EN 13770	Code	2568.804
Mounting device for socks necessary for Code 2568.804	Code	2568.806
Sample holder for ABRASION test on leather/coated fabrics according to UNI EN ISO 17076-2	Code	2568.260
ISO 9001 Calibration report	Code	2560.CC1
ISO 17025 Calibration Certificate (ACCREDIA)		

UPGRADING OPTIONAL ACCESSORIES ONLY FOR 6 POSITIONS

Optional additional position Ø 38 mm set, necessary to add an additional position to perform tests of ABRASION according to UNI EN ISO 12947-1 and PILLING according to ASTM D4970	Code	2568A.800
Optional additional position Ø 90 mm set, necessary to add an additional position to perform tests of PILLING according to UNI EN ISO 12945	Code	2568A.802
Set of common parts for each additional position, necessary to set-up any additional position, regardless from the type of sample holder matched	Code	2568A.818

REFERENCE STANDARDS

- ABRASION Ø 38 mm: UNI EN ISO 5470-2, UNI EN ISO 12947, UNI EN 13770, UNI EN ISO 20344, ISO 13520, ISO 17704, UNI EN ISO 17076-2, UNI EN 388 6.1.3 (2017), UNI EN 530, UNI EN 13520, ASTM D4966, BS 3424-24 TM 27A, UNI EN ISO 12947-1-2-3-4 (formerly BS 5690:1991), UNI EN ISO 12947-1-2-3-4 (formerly SN 198529), UNI EN ISO 12947 (formerly SFS 4328 - EMPA Method), IWTO 40, IWS TM 112, VDA 230-211, M&S P18C, M&S P19, M&S P19A, M&S P19B, M&S P19C, JIS L1018 6.18.5 meth E, ADIDAS GE 63.
- PILLING Ø 38 mm: IWS TM 196, ASTM D4970, M&S P140, M&S P17.
- PILLING Ø 90 mm: UNI EN ISO 12945-2, UNI EN ISO 12945-2 (formerly SN 198525).
- Officially approved by Marks & Spencer for:
ABRASION test: M&S P18C, M&S P19, M&S P19A, M&S P19B, M&S P19C;
PILLING test: M&SP17, M&S P140.

DIMENSIONS / POWER SUPPLY

Weight: 100 kg
Dimensions: (L) 800 x (W) 750 x (H) 410 mm
Power supply: 115 Vac or 230 Vac, 50/60 Hz, single-phase

Photographs and descriptions of the present leaflet have to be considered as purely indicative and not binding